

**MEDIMORE**  
**MEDICAL SUPPLIES**

# **DIAGNOSTIC CATALOGUE**

**Suppliers of Medical Disposables and Equipment**

**Direct Importers & Distributors of all medical equipment,  
consumables and surgicals.**

**We have Supplies. Do you have the *Patients*.**

## PREGNANCY TESTS

Our pregnancy tests are highly accurate and sensitive containing Colloidal Gold™ patented technology.



### FEATURES:

- World Health Standard
- High accuracy and sensitivity
- Excellent for retail shops and Pharmacies
- Used by Hospitals & Professionals worldwide.

AVAILABLE AT MOST WHOLESALERS AND PHARMACIES

PROD. CODE	DESCRIPTION	PACK SIZE	CASE SIZE	NAPPI CODE	BARCODE
CHPTS25	Pregnancy Test Strip CH - individually packed (R)	25	250	1004139001	6009701458275
CHPTS100	Pregnancy Test Strip CH - individually packed (R)	100	50	1004139001	6009701458282
CHPTD50	Pregnancy Device CH - individually packed (R)	50	36	1004137001	6009701458268
EPT001	Pregnancy Midstream Test - CH	1	400	184625001	0700371643843
EMPT25	Pregnancy Mainstream Test - CH (Shipper)	25	60	184625001	0700371643850
EPT002	Pregnancy Test - Double Pack - CH	1	500	184626001	0700371643836
LHOVA5	Ovulation Test - CH - (5 PACK)	5	20		0700371643867



## DRUG TESTS

Clinihhealth certified drug tests are rapid and indicates results in minutes.



### FEATURES:

- Highly sensitive drug test with results in minutes
- SAHMSA & FDA Approved
- Ideal for rehab centres, Doctors, Pharmacies and private use
- Tests for all major drug types including Cocaine, Dagga, Methamphetamine, Heroin, TIK, Benzodiazepine
- Methadone MDMA Ketamine Opiates

PROD. CODE	DESCRIPTION	PACK SIZE	CASE SIZE	NAPPI CODE	BARCODE
D5CH1	Drug Test 5 Clinihhealth (Single)	1	500	210491001	6009879738292
D6CH1	Drug Test 6 Clinihhealth (Single)	1	500	210566001	6009879738308
D6CH20	Drug Test 6 Clinihhealth 20s	20	10		6009701450972
D8CH20	Drug Test 8 Clinihhealth 20s	20	10		6009701450965
DRUGCUP6	Drug Test Clinihhealth 6 panel CUP	1	50	274344001	6009701452280
DRUGCUP12	Drug Test Clinihhealth 12 panel CUP	1	50		6009701451009
MQL20	Methaqualone Test 25s	1	1		
CHTHC50	Dagga/THC testClinihhealth 50s	1	10	274342001	6009701450705
SALDT6CH	Saliva Drug Test 6	1	20		6009701452648



## INFECTIOUS DISEASE TESTS

Clinihealth infectious disease tests are high quality certified tests trusted globally

- Quick Results
- Swabs & Lancets Included
- High Sensitivity



AVAILABLE  
IN PACK  
of 25s

### FEATURES:

- Rapid point of care
- High quality products
- Manufactured by Certified Laboratories
- High Specificity

PROD. CODE	DESCRIPTION	PACK SIZE	CASE SIZE	NAPPI CODE	BARCODE
CHSYPC25	Syphilis Test 25s	25	50		6009701450736
CHHIV1	HIV Kit 1 & 2 Compl. Clinihealth	1	600	27435001	6009879738377
CHHIV25	HIV Kit Bulk - 25s	25	50		6009701452785
CHMAL25	Malaria Test Device Clinihealth	25	10	196710001	6009701451870
CHMAL1	Malaria Test Device Complete Clinihealth	1	500	198398001	6009879738385



## URINE ANALYSIS

Quality care with trusted products from Clinihealth®

### FEATURES:

- Highly sensitive metabolic tests for professional use
- Each bottle contains 100 strips with manufacture seal of quality for point of care
- Analyses Leukocytes / Nitrites / Glucose / Protein / PH / Blood / Uribilogen / Specific Gravity / Ascorbic acid / Ketone / Bilirubin.



### Urine Analyser Features

Quick Automation  
High Quality Analyser used for Quick and Rapid Reading of Combi Strips  
Test Principle : Super High luminosity cold light source reflection  
Test Parameters : up to 11  
Test Speed Up to 120 Tests per hour on continuous Mode  
Memory up to 1000  
Absorbs Excess Urine Automatically - Prevents Cross infection  
Thermal Printer inside with outer interface  
Connects to pc - 240V - Weight 4kg  
Dimensions : 355mmx300mmx145mm  
**Warranty Included**



PROD. CODE	DESCRIPTION	PACK SIZE	CASE SIZE	NAPPI CODE	BARCODE
URICH6P	Urinalysis 6P Clinihealth	100	100	194717001	6009701458299
URI10A	Urinalysis Combi 10 Clinihealth	100	100	274359001	6009879738407
URIANL1	Urine analyser	1	1		

## FIRST AID KITS

First Aid Kits in South Africa are regulated by the occupational health and safety act which states that an employer shall take necessary steps to ensure employees receive First Aid treatment from a trained first aider in case of accidents and emergencies.

First Aid is an important skill. By performing simple procedures and following certain guidelines, it may be possible to save lives by giving basic treatment until professional medical help arrives.



PROD. CODE	DESCRIPTION	PACK SIZE	BARCODE
REG3REFW	Regulation 3 Refill wooden splint	1	6009701455373
FAREG3M	Complete in metal box Wall Mountable	1	6009701455380
FARREG3P	Complete in Plastic utility box First Aid Reg 7	1	6009701455397
REG7REFW	Regulation 7 Refill Wooden Splint	1	6009701455434
FAREG7M	Complete in Metal Box - wall mountable	1	6009701455441
FARREG7P	Complete in Plastic Utility Box	1	6009701455458
FACLI303	First Aid Kit - 303 pieces in zip Bag - Clinihealth	1	6009701450767
FACLI210	First Aid Kit - 210 pieces in zip Bag - Clinihealth	1	6009701450767
PULSEOXIPC60	Pulse Oximeter PC60	1	6009701450163
REFILLREG3	Refill Regulation 3	1	6009701450187
REFILLREG7	Refill Regulation 7	1	6009701450194



## DISPOSABLE WEAR RANGE

Leading the way to better healthcare

### FEATURES:

- High quality graded material
- Easy comfort and care
- Various colours and sizes
- We can custom make any apron type, micron, colour and pack sizes.



LINEN SAVERS




MOB CAPS



APRONS



PROD. CODE	DESCRIPTION	PACK SIZE	BARCODE
APRPE100ECO	Aprons Econo Disposable - WHITE	100	6009701455014
APR20MC100WHT	Aprons 20mc Disposable - WHITE	100	6009701455021
APR20MC100CLR	Apron 20mc Disposable CLEAR/COLOURS	100	6009701455038
MOBCAP100	Mob Cap - Single Elastic 100s	100	6009701455045
LIN2006PLY	Linen Savers 200 - 6 Ply	100	6009701453069
LIN2004PLY	Linen Savers 200 - 4 Ply	100	6009701455076
SMK20MC100	Apron Smock with Sleeve - various colours	100	6009701450200



MEDIMORE  
MEDICAL SUPPLIES